


Main Prosthetic Components



Healing

Closure Screw



Healing Cap - Emergency Ø3.2



Healing Cap - Emergency Ø3.6




Healing Cap - Emergency Ø5.0




Impression


Impression Coping



Analog



Main Prosthetic Components



Screw Retained Restoration

Castable



Cobalt-Chrome Castable



Temporary Abutment



Angled Cobalt-Chrome



ANTI ROT

MiP 3.0 NP 3.5

VCON04S0001 VCON04N0001

VCON04S0002 VCON04N0002

ROT

MiP 3.0 NP 3.5

VCON05S0001 VCON05N0001

VCON05S0002 VCON05N0002

MiP 3.0 NP 3.5

VCON06S0001 VCON06N0001

VCON06S0002 VCON06N0002

17° 30°

1 VCHA170001 1 VCHA300001

2 VCON05N0003 2 VCON05N0003

17° 30°

1 VCHA170001 1 VCHA300001

2 VCON05N0004 2 VCON05N0004

Cemented Prosthesis



Clinical Screw



h1.5mm: VCON07N1501

h2.5mm: VCON07N2501

h3.5mm: VCON07N3501

h4.5mm: VCON07N4501

15° 25°

15° h1.5 mm: VCON07N1515 25° h1.5 mm: VCON07N1525

15° h2.5 mm: VCON07N2515 25° h2.5 mm: VCON07N2525

15° h3.5 mm: VCON07N35015 25° h3.5 mm: VCON07N3525

MiP 3.0 NP 3.5

TITANIUM 5 TITANIUM 5

VCON09S07 VCON09N07

max 15 Ncm max 25 Ncm

TITANIUM DLC TITANIUM DLC

VCON09S07T VCON09N07T

Angled Prosthetic Screws

NP 3.5

TITANIUM 5

VCON09N00-TLB

max 25 Ncm

Short VSDTLB-1

Medium VSDTLB-2

Long VSDTLB-3